

# RK-05

**Artikelnummer:** RK-05



## Short Description

Portable rack

Tilt and stay feature (30 degrees)

Fast assembly

## Beschreibung

If you want to observe your object or specimen from a different perspective, the Dino-Lite RK-05 stand provides a multi-angle view with its 30 degrees inclinable pole. RK-05 is a light-weight, ultra-portable stand with a one-step assemble/disassemble design.

All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-05F with fine focus adjustment](#)

## Spezifikation

| Compatibility    |                                                           |
|------------------|-----------------------------------------------------------|
| Compatible with: | All Dino-Lite models, except Long Working Distance models |
| Housing          |                                                           |

|                      |                                  |
|----------------------|----------------------------------|
| Dimensions           | 14cm (W) x 11cm (L) x 16.5cm (H) |
| Weight               | 611 g                            |
| Features             |                                  |
| ESD safe             | No                               |
| Information          |                                  |
| Package contents     | Stand, HD-P1 holder              |
| Warranty information | 2 years European warranty        |
| Accessory range:     | Professional stands              |
| Price range          | €150,00 - €250,00                |

## Product Gallery

