

# RK-05F

**Artikelnummer:** RK-05F



## Short Description

Portable rack  
Tilt and stay feature (30 degrees)  
Fast assembly  
Fine-focus adjustment

## Beschreibung

RK-05F is a light-weight, ultra-portable stand with a one-step assemble/disassemble design.

If you want to observe your object or specimen from a different perspective, the Dino-Lite RK-05F stand provides a multi-angle view with its 30 degrees inclinable pole feature. Additionally it has an integrated fine focus adjustment knob to make focusing easier and more precise.

RK-05F is a light-weight, ultra-portable stand with a one-step assemble/disassemble design.

All Dino-Lite microscopes are compatible with this stand with the exception of the Long Working Distance microscopes.

The package includes the stand and a HD-P1 holder.

[Also available: RK-05 without fine focus adjustment](#)

## Spezifikation

| Compatibility        |                                                           |
|----------------------|-----------------------------------------------------------|
| Compatible with:     | All Dino-Lite models, except Long Working Distance models |
| Housing              |                                                           |
| Dimensions           | 14cm (W) x 11cm (L) x 16.5cm (H)                          |
| Weight               | 636 g                                                     |
| Features             |                                                           |
| ESD safe             | No                                                        |
| Information          |                                                           |
| Package contents     | Stand, HD-P1 holder                                       |
| Warranty information | 2 years European warranty                                 |
| Accessory range:     | Professional stands                                       |
| Price range          | €200,00 - €300,00                                         |

## Product Gallery

