

# CA1070

**Product Code:** CA1070



## Short Description

Sturdy watertight carrying case

## Description

Sturdy, watertight case for Dino-Lite microscopes. This heavy duty carrying case with foam inlay protects your Dino-Lite against dust, scratches and moisture. Watertight according to IP67 standard.

## Specification

| Compatibility        |                                                       |
|----------------------|-------------------------------------------------------|
| Compatible with:     | All Dino-Lite models                                  |
| Housing              |                                                       |
| Dimensions           | Exterior: 210 x 167 x 90mm Interior: 186 x 123 x 75mm |
| Weight               | 570 g                                                 |
| Features             |                                                       |
| ESD safe             | No                                                    |
| Information          |                                                       |
| Package contents     | Case, foam inlay                                      |
| Warranty information | 2 years European warranty                             |
| Accessory range:     | Various                                               |
| Price range          | €0 - €50,00                                           |

Product Gallery

