

DermaScope® Wireless (MEDLW4DW)



Código del producto: MEDLW4DW

Short Description

1.3 Megapixel resolution
Versatility due to exchangeable caps
10-50x Magnification
Wireless
Medical Device Class 1

Descripción

The Dino-Lite DermaScope® Wireless (MEDLW4DW) has a 1.3 megapixel camera that can capture sharp images of the skin and features a built-in and fully adjustable polarization filter to greatly reduce the gloss effect of the skin.

The Dino-Lite DermaScope® Wireless is suitable for any physician or dermatologist or skin specialist. For those who put even more demands on the equipment and want to have the latest features and capabilities, the DermaScope® Polarizer HR is also available (MEDL7DW)




USB 2.0



Exchangeable
cap



Adjustable
~10-50x



Especificación

Lighting	
Light/ LED type	White
Number of LEDs	8
LED on/off switchable:	Yes
Infrared filter	IR cut-filter >650 nm
Diffuser available	No
Polarizer	Yes, linear
Optics	
Magnification	10-50x
Lens type	Glass with anti-reflection coating
Sensor	
Sensor type	CMOS
Resolution	1.3 Megapixel (1280x1024)
Maximum frame rate	30 fps
Compatibility	
Interface	USB 2.0
Operating system	Windows 7, 8, 10 & 11, MacOS 10.9 and up
Software	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
Supported image formats (Windows)	BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
Supported video formats (Windows)	WMV, FLV, SWF
Supported image formats (MacOS)	JPEG, PNG
Supported video formats (MacOS)	MOV
Imaging standards	DirectShow, UVC
Wifi	Yes
Housing	
Housing material	Composite/plastic housing
Magnification lock	Yes
Dimensions	10.5 cm (L) x 3.2 cm (D)

Weight	105 g
Cable length	1.8m
Features	
Special feature	No
Measurement	No
Calibration	No
Microtouch sensor	Yes
Information	
Package contents	Microscope, carry pouch, software CD, user manual, N3C-C, N3C-L
Warranty information	2 years European warranty
Regulatory approval	Medical Device Class 1 – Medical Devices Regulation (EU) 2017/745
Price range	€900,00 - €1100,00

Product Gallery

