

EarScope Basic Wireless (MEDLW3E)

Merk:

Productcode: MEDLW3E



Korte omschrijving

640x480 Resolution

Basic model for otoscopy

20-30x Magnification

Wireless

Medical Device Class 1

Omschrijving

The EarScope Basic Wireless (MEDLW3E) is a compact and safe solution to look at the ear. The cap is detachable for hygienical reasons. The MEDLW3E has a 640x480 resolution and magnification of 20 ~ 30x. This model is delivered with disposable specula of 3, 4 and 5 mm (10 each).



USB 2.0



Adjustable
~20x-30x



VGA resolution
640 x 480



4 White LEDs
non switchable



Wireless

Specificatie

Verlichting

Licht / LED-type

White

Aantal leds	4
LED aan / uit schakelbaar:	Yes
Infrarood filter	IR cut-filter >650 nm
Diffusor beschikbaar	No
Polarisator	No
Optiek	
Magnification	20-30x
Type lens	Glass with anti-reflection coating
Sensor	
Sensortype	CMOS
Resolutie	VGA (640x480)
Maximale framesnelheid	30 fps
Compatibiliteit	
Interface	USB 2.0
Besturingssysteem	Windows 7, 8, 10 & 11, MacOS 10.9 and up
Software	DinoCapture 2.0 (Windows) cable + wireless, DinoXcope (Mac OS) cable only, DinoConnect wireless only
Ondersteunde beeldformaten (Windows)	BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
Ondersteunde videoformaten (Windows)	WMV, FLV, SWF
Ondersteunde beeldformaten (MacOS)	JPEG, PNG
Ondersteunde videoformaten (MacOS)	MOV
Beeldvormingsnormen	DirectShow, UVC
Wifi	Yes
Behuizing	
Materiaal van de behuizing	Composite/plastic housing
Vergrotingsslot	No
Afmetingen	10.5cm (L) x 3.2cm (D)
Gewicht	95 g
Kabellengte	1.8m
Kenmerken	
Meting	No
Kalibratie	No
Microtouch-sensor	No
Informatie	
Inhoud van de verpakking	Microscope, carry pouch, 10x 3 mm specula, 10x 4 mm specula, 10x 5 mm specula, software CD, user manual

Garantie-informatie	2 years European warranty
Regelgevende goedkeuring	Medical Device Class 1 – Medical Devices Regulation (EU) 2017/745
Prijsklasse	€500,00 - €700,00

Product galerij

