

EarScope Pro Wireless (MEDLW4E)

Productcode: MEDLW4E



Korte omschrijving

1.3 Megapixel resolution

Pro model for otoscopy

55-90x Magnification

Wireless

Medical Device Class 1

Omschrijving

The EarScope Pro Wireless (MEDLW4E) is a compact and safe solution for watching the eardrum and the external ear. The MEDL4E is the Pro model with 1,3 Megapixel resolution and a magnification of 55 ~ 90x. This model is delivered with disposable specula of 3, 4 and 5 mm (10 each).



USB 2.0



6 White LEDs
switchable



Adjustable
~55x-90x



1.3 megapixel
1280 x 1024



Magnification
lock



Wireless

Specificatie

Verlichting

Licht / LED-type	White
Aantal leds	6
LED aan / uit schakelbaar:	Yes
Infrarood filter	IR cut-filter >650 nm
Diffusor beschikbaar	No
Polarisator	No
Optiek	
Magnification	55-90x
Type lens	Glass with anti-reflection coating
Sensor	
Sensortype	CMOS
Resolutie	1.3 Megapixel (1280x1024)
Maximale framesnelheid	30 fps
Compatibiliteit	
Interface	USB 2.0
Besturingssysteem	Windows 7, 8, 10 & 11, MacOS 10.9 and up
Software	DinoCapture 2.0 (Windows) cable + wireless, DinoXcope (Mac OS) cable only, DinoConnect wireless only
Ondersteunde beeldformaten (Windows)	BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
Ondersteunde videoformaten (Windows)	WMV, FLV, SWF
Ondersteunde beeldformaten (MacOS)	JPEG, PNG
Ondersteunde videoformaten (MacOS)	MOV
Beeldvormingsnormen	DirectShow, UVC
Wifi	Yes
Behuizing	
Materiaal van de behuizing	Composite/plastic housing
Vergrotingsslot	Yes
Afmetingen	10.5cm (L) x 3.2cm (D)
Gewicht	95 g
Kabellengte	1.8m
Kenmerken	
Bijzonderheid	No
Meting	No
Kalibratie	No
Microtouch-sensor	Yes
Informatie	
Inhoud van de verpakking	Microscope, carry pouch, 10x 3 mm specula, 10x 4

	mm specula, 10x 5 mm specula, software CD, user manual
Garantie-informatie	2 years European warranty
Regelgevende goedkeuring	Medical Device Class 1 – Medical Devices Regulation (EU) 2017/745
Prijsklasse	€800,00 - €1000,00

Product galerij

