

AM4025X

Merk: Dino-Lite

Productcode: AM4025X



Korte omschrijving

1.3 Megapixel resolution

Fits 23, 30 and 30,5 mm ocular

Advanced software

Including measurements

Omschrijving

The DinoEye Edge 1.3MP eyepiece camera provides superior image quality and true colour reproduction. This makes this DinoEye suitable for demanding applications where true colour reproduction is of the utmost importance, like with pathology samples or quality control.



USB 2.0



Measurement
functionality



1,3 Megapixels
1280 x 1024



connects
to c-mount



Fits 23, 30,
30,5 mm oculars



Extended
Depth Of Field
EDOF

The large field of view allows you to capture more details with one image. You can replace one of the eyepieces of the microscope by the Dino-Eye camera which has a diameter of 23mm, or you can use one of the provided adapters of 30mm or 30,5mm to make it fit.

Additionally you can use the c-mount thread that is included with the 30,5mm adapter to be able to connect the DinoEye to the c-mount adapter (not included) on the Trinocular tube of the microscope. The AM4025X eyepiece camera will therefore fit to almost all microscope brands and models available. The AM4025X comes with DinoCapture software for Windows and DinoXcope software for MAC, which includes measurement, annotation, video and time-lapse functions.

Specificatie

Optiek	
Magnification	20x, built-in
Type lens	Glass with Anti-Reflection coating
Sensor	
Sensortype	CMOS
Resolutie	1.3 Megapixel (1280x1024)
Maximale framesnelheid	30 fps
Compatibiliteit	
Interface	USB 2.0
Besturingssysteem	Windows 7, 8, 10 & 11, MacOS 10.12 and up
Software	DinoCapture 2.0 (Windows), DinoXcope (Mac OS)
Ondersteunde beeldformaten (Windows)	BMP, GIF, PNG, JPG, TIF, RAS, PNM, TGA, PCX, MNG, WBMP, JP2, JPC, PGX
Ondersteunde videoformaten (Windows)	WMV, FLV, SWF
Ondersteunde beeldformaten (MacOS)	JPEG, PNG
Ondersteunde videoformaten (MacOS)	MOV
Beeldvormingsnormen	DirectShow, UVC
Wifi	Wireless-ready, requires the WF-10 WiFi streamer (optional)
Behuizing	
Materiaal van de behuizing	Metal body
Afmetingen	92 mm (H) x 32 mm (D)
Gewicht	80 g
Kabellengte	1.8 m
Ocular outer diameter:	23, 30, 30,5 mm
Mounting type:	Inside ocular/tube & c-mount
Kenmerken	
Bijzonderheid	Extended Depth of Field (EDOF)
Meting	Yes

Kalibratie	Yes
Microtouch-sensor	No
Informatie	
Inhoud van de verpakking	DinoEye eyepiece camera, carry pouch, software CD, calibration target, user manual
Garantie-informatie	2 years European warranty
Regelgevende goedkeuring	CE, FCC, ROHS
accessoires:	DinoEye Eyepieces
Prijsklasse	€250,00 - €500,00

Product galerij

