

# RK-02

Ürün Kodu: RK-02



## Short Description

Gooseneck stand

## Açıklama

The RK-02 is a durable gooseneck stand solution designed for flexible and versatile usage with a Dino-Lite digital microscope. The C-clamp of the RK-02 fits to virtually any surface of up to 7 cm (~3 inch) thick, and can be mounted horizontally or vertically. With this sturdy stand, you can easily position the Dino-Lite by bending the 60 cm long arm and adjusting the rotating head. The base of the RK-02 can also be rotated providing even more freedom and it can be fixed with the quick release fastener for use on uneven surfaces. The RK-02 gooseneck arm is the ideal choice to flexibly use any of the Dino-lite models, except the high magnification models.

## Özellikler

| Compatibility    |                                                        |
|------------------|--------------------------------------------------------|
| Compatible with: | All Dino-Lite models, except High Magnification models |
| Housing          |                                                        |
| Dimensions       | Gooseneck 60cm (L) x 1.5cm (D)                         |
| Weight           | 1.5 kg                                                 |
| Features         |                                                        |
| ESD safe         | No                                                     |

| Information          |                           |
|----------------------|---------------------------|
| Package contents     | Stand, HD-P1 holder       |
| Warranty information | 2 years European warranty |
| Accessory range:     | Professional stands       |
| Price range          | €100,00 - €200,00         |

## Product Gallery

